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<https://www.molnlycke.com/zflo/>

Adult

Z-Flo™ is a mark owned by EdiZONE, LLC of Alpine, Utah

Mölnlycke® Z-Flo™ Fluidized positioner

Mold the positioner to desired shape

The positioner maintains its shape

en Please read all of the instructions before use. Use as directed by appropriate clinician. Always follow your facility protocol for turning and positioning.

Mölnlycke® Z-Flo™ Fluidized positioner Adult

Positioner product for Adult and Adolescent hospital population.

Product Description

The positioner consists of a fluidized medium composed of silicone oil, microspheres and polyethylene beads which is encapsulated in a polyurethane film exterior. The positioner can be molded to support multiple anatomical shapes and clinical needs.

Intended use

Mölnlycke® Z-Flo™ Fluidized positioner is designed to provide patient positioning, pressure redistribution and/or to offload anatomical sites, which may be used as part of a prophylactic therapy to help prevent pressure ulcers.

Instructions for use

Adult and Adolescent

- Patient positioning, turning and repositioning
- Off-loading of heel, sacrum, occiput and other bony prominences
- Contracture management

Lifetime

Product properties will deteriorate over time in use. The rate of deterioration is dependent on use environment and application. The product can be used with one patient until the product has deteriorated and does no longer fulfil its intended use.

Other information

If the product is used after expiry date product properties cannot be ensured. The product expiry date is printed on the packaging labels.

Instructions for use

1. Remove the positioner from the transparent plastic bag.
2. Positioner cover or liner is required between the patient's skin and the positioner at all times during use.
3. Place the positioner appropriately under the patient.
4. Gently mold the positioner towards the patient/anatomical site to support and maintain the desired position, achieving pressure redistribution to minimize excessive forces on the skin.
5. The positioner can be molded to accommodate medical equipment or airway clearance.
6. The positioner can easily be remolded to accommodate changed clinical needs.
7. If the product has deteriorated (e.g. due to time in use) and no longer has the ability to hold the desired position it should be discarded and replaced.

The positioners are MR Safe and can also be used during CT. If used in connection to X-ray, do not X-ray through the positioner as it may distort the image.

Manufacturer

Mölnlycke Health Care AB
Gamlestadvägen 3C, Box 13080, SE-402 52 Göteborg, Sweden

Patent numbers:
US 8,171,585 B2
US 9,120,666 B2

Toll free number:
USA 1-800-882-4582
Canada 1-800-494-5134

Australian Sponsor Address:
Mölnlycke Health Care Pty. Ltd.
12 Narabang Way, Belrose, NSW 2085, Australia

MD Medical Device

40°C

104°F

Keep dry
Conserver au sec

Issued 2021-04

Master PD-54595 rev 09
PD-540085 rev 04

Made in USA

The Mölnlycke trademark, name and logo are registered globally to one or more of the Mölnlycke Health Care Group of Companies

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